



Protocol: IV-HSL Emitter Module

Printable Version

Overview

This protocol reconstitutes the BjaI/BjaR quorum sensing components from *Bradyrhizobium japonicum* to establish IV-HSL-producing synthetic cells (emitters) and IV-HSL-responsive *Escherichia coli* cells (receivers), implementing the [IV-HSL Emitter Cell](#).

BjaI is expressed inside Emitter Cells containing PURExpress to produce the enzyme BjaI from the template `pT7-bjaI`. BjaI will catalyze a reaction between the membrane impermeable IV-CoA and SAM substrates to yield membrane permeable IV-HSL.

E. coli cells expressing BjaR act as receiver cells, providing an easy means to detect IV-HSL production. When BjaR binds IV-HSL, expression of a fluorescent reporter gene controlled by a BjaR-regulated promoter is triggered.

Successfully built IV-HSL Emitter Cells will release IV-HSL and induce GFP expression in XL10-Gold cell with increasing green fluorescence over time.

There are five key stages to making the IV-HSL Emitter Cell:

| Step | Process | Hands-on Time | Total Time | Notes |
|------|--|---------------|------------|--|
| 1 | <u>Pre-culture BjaR receiver cells</u> | 30 mins | 3.5 hr | |
| 2 | <u>Prepare lipids-in-oil solution, outer solution, and substrate stock solutions</u> | 1 hr | 4 h | Buffers and lipids may be prepared in advance and used for experiments on subsequent days. |
| 3 | <u>Assemble PURE reactions</u> | 30 mins | 30 mins | |
| 4 | <u>Encapsulate liposomes</u> | 30 mins | 30 mins | |
| 5 | <u>Measure and image</u> | 30 mins | 6–12 h | Total time depends on the exact experiment and incubation conditions. GFP expression should be seen over the first 6 hours at 37C. |

Materials and Equipment

| Name | Product | Manufacturer | Part # | Price | Link |
|-----------------------------|--|-------------------|------------------|----------|------------------------|
| Buffers | | | | | |
| Glucose | D-(+)-Glucose, 99% | Thermo Scientific | A16828-36 | \$41.65 | [link] |
| Sucrose | Sucrose, 99% | Thermo Scientific | A15583-36 | \$41.65 | [link] |
| Lipids | | | | | |
| Egg PC | 25mg/mL | Avanti Lipids | 840051C-200mg | \$186 | [link] |
| Liss-Rhod-PE | 18:0 Liss Rhod PE 1 mg/mL | Avanti Lipids | 810179P-1mg | \$273.47 | [link] |
| Mineral Oil | Mineral oil, mixed weight | Thermo Scientific | AC415080010 | \$53.40 | [link] |
| Glass Syringe 250 uL | | Hamilton | 14-815-238 | \$150.15 | [link] |
| PURE | | | | | |
| PURE | PURExpress | NEB | E6800S | \$295.00 | [link] |
| RNase Inhibitor | RNase Inhibitor, Murine | NEB | M0314S | \$81.00 | [link] |
| DNA | pT7-bjaI | b. next | | | [link] |
| | bjaR-GFP-native | b.next | | | [link] |
| OptiPrep | OptiPrep - Density Gradient Media (Iodixanol) | COSMO BIO USA | AXS-1114542 | \$172 | [link] |
| SAM | S-adenosylmethionine (SAM) | NEB | B9003S | \$45 | [link] |
| IV-CoA | Isovaleryl coenzyme A lithium salt hydrate | Millipore Sigma | I9381-10MG | \$348 | [link] |
| IV-HSL | 3-Methyl-N-[(3S)-tetrahydro-2-oxo-3-furanyl]butanamide | LGC | TRC-M282980-50MG | \$171 | [link] |

| Name | Product | Manufacturer | Part # | Price | Link |
|------------------------|--|-------------------|-----------|-------|------------------------|
| DMSO | Dimethyl sulfoxide | Thermo Scientific | 042780.M1 | \$342 | [link] |
| | | | | | |
| Cell culture | | | | | |
| XL10-Gold Cells | XL10-Gold Ultracompetent Cells | Agilent | 200314 | \$223 | [link] |
| M9 Media | M9, Minimal Salts, 5X, powder, minimal microbial growth medium | Sigma-Aldrich | M6030-1KG | \$260 | [link] |

Step 1: Pre-culture BjaR receiver cells

- ☐ Prepare glycerol stock of BjaR receiver cells
 - ☐ Transform XL-10 Gold competent *E. coli* with **bjaR-GFP-native** :
 - ☐ Add 1–5 µl containing 1 pg–100 ng of plasmid DNA **bjaR-GFP-native** to 50 µl of XL10-Gold cell mixture. Carefully flick the tube 4–5 times to mix cells and DNA. **Do not vortex.**
 - ☐ Place the mixture on ice for 15 minutes. Do not mix.
 - ☐ Heat shock at exactly 42°C for 40 seconds. Do not mix.
 - ☐ Place on ice for 5 minutes. Do not mix.
 - ☐ Pipette 950 µl of room temperature SOC into the cell mixture.
 - ☐ Shake the cell mixture vigorously (250 rpm) at 37°C for 60 minutes.
 - ☐ Warm Ampicillin LB agarose plates at 37°C for 10 mins.
 - ☐ Mix the cells thoroughly by flicking the tube and inverting, then perform several 10-fold serial dilutions in LB.
 - ☐ Spread 50–100 µl of each dilution onto a Ampicillin agarose plate and incubate overnight for ~15 hrs at 37°C.
 - ☐ [if we're including making a glycerol stock, need overnight culture and glycerol stock preparation here]
- ☐ Prepare a streak plate from the glycerol stock ([reference](#))
 - ☐ Streak a Ampicillin LB plate from the glycerol stock and incubate overnight at 37C.
- ☐ Prepare M9 Media containing 1× M9 salts, 0.34 mg/ml–1 thiamine hydrochloride, 0.2% casamino acids, 2 mM MgSO₄, 100 µM CaCl₂ and 0.4% (wt/vol) glucose.

- ☐ Pick a colony from the *E. coli* streak plate, and inoculate a 5 mL culture tube containing the M9 media with 100 ug/mL carbenicillin.
- ☐ Incubate the cells at 37 °C, 225 rpm, for 3 h. *Prepare Emitter liposomes while the cells incubate.*
- ☐ Dilute the culture media with the pre-warmed M9 media until OD600 = ~0.1.
- ☐ Balance osmolarity of the culture media with PURE (inner solution in liposomes) by adding glucose to the M9 media:

| | Volume to mix (uL) |
|-------------------|--------------------|
| M9 media | 1000 |
| 3M Glucose | 293.81 |

Step 2: Prepare lipids-in-oil solution, outer solution, and substrate stock solutions

Prepare lipids-in-oil (mineral oil) solution

- ☐ Clean glass syringes.
 - ☐ Pour a small amount of 95% ethanol into a glass container (e.g. a 10 mL beaker).
 - ☐ Assemble the glass syringe and prime it by drawing ethanol into the glass syringe, then empty into a waste bottle.
- ☐ Use glass syringes to add lipids, as shown in the table below, into the 10 ml glass vial containing 1 ml of mineral oil (final lipid concentration is 5 mg/ml).

| Lipids | Stock Concentration (mg/mL) | Volume to add (uL) | Target percentage |
|--------------------------|-----------------------------|--------------------|-------------------|
| Egg PC | 25 | 160 | 66.68 |
| Cholesterol | 50 | 20 | 33.32 |
| 18:0 Liss Rhod PE | 1 | 5 | 0.01 |

- ☐ Heat the lipids-in-oil mixture on a hotplate at 55 C for 3 hrs.
- ☐ Vortex the lipids-in-oil mixture for 1 min.
 - The lipids-in-oil mixture can be stored at 4 C for up to 3 days.

Prepare outer solution

Final concentration of sugar stock solution is 900 mM


| Buffer | Volume to add (uL) |
|-------------------------|--------------------|
| 3M Glucose Stock | 700 |
| H2O | 300 |

Prepare substrate stock solutions

| Substrate | Concentration (uM) | MW (g/mol) | Weight (g) | Final Volume (mL) |
|---------------|--------------------|------------|------------|-------------------|
| SAM | 5000 | 398.44 | 1.99 | 1 |
| IV-CoA | 5000 | 851.65 | 4.26 | 1 |
| IV-HSL | 10 | 183.21 | 1.83 | 1 |

Step 3: Assemble PURE Reactions


PURE reaction setup

| | Sample | Negative control | Positive control | |
|---|-------------|------------------|------------------|---|
| Component | Volume (uL) | Volume (uL) | Volume (uL) | Notes |
| PURE Solution A | 12 | 12 | 0 | PURE energy solution: small molecules |
| PURE Solution B | 9 | 9 | 0 | PURE proteins and ribosomes |
| RNAse Inhibitor | 1.5 | 1.5 | 0 | Prevents RNAse activity |
|  <u>EM01-pOpen-pT7-BjaI</u> (~200 ng/uL) | 1.5 | 0 | 0 | DNA encoding green fluorescent protein |
| SAM (5mM) | 1.8 | 1.8 | 0 | Substrate for IV-HSL production. |
| IV-CoA (5mM) | 0.48 | 0.48 | 0 | Substrate for IV-HSL production. |
| OptiPrep | 1.5 | 1.5 | 1.5 | Adds density for phase-transfer |
| IV-HSL (10 uM) | 0 | 0 | 0.3 | Commercial IV-HSL for positive control. |
| 3M Glucose | 0 | 0 | 8.46 | |
| ddH2O | 2.22 | 3.72 | 19.74 | |
| Total | 30 | 30 | 30 | |

☐ Thaw reagents on ice and then keep on ice.

☐ Prepare a PCR strip in a strip holder on ice for assembly of the three reactions (Sample, Negative, Positive).

Step 4: Encapsulate PURE reactions into Liposomes

Some tips and tricks can be found in  ["Hello, world" PURE Liposomes](#).

- ☐ Set up a microfuge tube rack, with three 1.5 mL microfuge tubes per liposome encapsulation:
 - ☐ Number the tubes per the number of reactions assembled in Step 3.
 - ☐ For each reaction, label the two tubes:
 - ☐ **I** — Oil emulsion
 - ☐ **O** — Outer solution
- ☐ Add 30 ul of PURE reactions prepared in **Step 3** to tubes labelled **I**.
- ☐ Add 180 uL of the lipids-in-oil mixture on top of the PURE reactions in tubes labelled **I** and pipette vigorously until the emulsion becomes cloudy.
- ☐ Add 300 uL of outer solution to each of the tubes labelled **O**.
- ☐ Add 210 uL of the milky solution carefully on top of the outer solution in the tubes labelled **O**.
- ☐ Centrifuge at 9000 rpm at 4c for 10 mins.
- ☐ Remove the top oil and resuspend the pellet in 100 ul of outer solution.
- ☐ Collect the liposomes.

Step 5: Measure and Image Liposomes and Cells

Imaging using confocal microscopy (Operetta CLS):

While microscopy setups may vary, our performance data was collected using the following configuration.

- ☐ Add BjaR receiver cells prepared in Step 1 into 384 Well Glass Bottom Microplates.
- ☐ Add 10 uL of liposomes made in Step 3 on top of the receiver cells in 384 Well Glass Bottom Microplates.
- ☐ Imaging conditions using Operetta:
 - Temperature: 37 C degree
 - Green fluorescence channel (200 us exposure 95%) - excitation: 460-490 nm; emission: 500-550 nm.
 - Red fluorescence channel (50 us exposure 95%) - excitation: 530-560 nm; emission: 570-650 nm.
 - Brightfield (20 us 95%)
 - We capture a 6 h time lapse with 10 min intervals.
 - We also acquired z-stack images spanning from 0 μm to 80 μm of the focal plane.

Measuring using plate reader (BioTek Cytation 5):

- ☐ Add BjaR receiver cells prepared in Step 1 into 96 Well Glass Bottom Microplates.
- ☐ Add 10 uL of liposomes made in Step 3 on top of the receiver cells in 96 Well Glass Bottom Microplates.
- ☐ Procedures:
 - Temperature: 37 C degree
 - Read the fluorescence intensity from the bottom
 - Excitation wavelength: 485 nm ; Emission wavelength: 528 nm
 - We capture a 6 h time lapse with 5 min intervals

Background Protocols

- ☐ Prepare lipids for use in encapsulation: ➡ [Lipid Preparation](#)
- ☐ Prepare inner and outer buffers: 💧 [PURE inner and outer solution](#)

Resources and References

• Other Protocols

- Transformation protocol: 🧫 [Transformation v1.0.](#)

• Papers

- Smith, J. M., Hartmann, D. & Booth, M. J. Engineering cellular communication between light-activated synthetic cells and bacteria. *Nature Chemical Biology* **19**, 1138–1146 (2023).
[\[https://www.nature.com/articles/s41589-023-01374-7\]](https://www.nature.com/articles/s41589-023-01374-7)

Credits

- b.next